

[All Products](#) / [Lab and Production](#) / [Optical Cleaning](#) / [Gloves, Brushes, and Swabs](#)

Cotton-Tipped Swab Applicators



- Lint-Free
- 3" or 6" Wooden Stems
- Bags, Boxes, and Cases Available

Technical Information

Products

Title	Stock Number	Price	Buy
3" Length Bag, Cotton-Tipped Swabs	#56-926	A\$9.20	14 In Stock
3" Length Box (10 Bags), Cotton-Tipped Swabs	#56-927	A\$43.20	14 In Stock
3" Length Case (10 Boxes), Cotton-Tipped Swabs	#56-928	A\$240.00	Contact Us
6" Length Bag, Cotton-Tipped Swabs	#56-924	A\$9.20	20+ In Stock
6" Length Box (10 Bags), Cotton-Tipped Swabs	#39-944	A\$43.20	20+ In Stock
6" Length Case (10 Boxes), Cotton-Tipped Swabs	#56-925	A\$240.00	5 In Stock



Copyright 2024, Edmund Optics Singapore Pte. Ltd, 18 Woodlands Loop #04-00, Singapore 738100

[PRIVACY POLICY](#) [COOKIE POLICY](#) [TERMS & CONDITIONS](#)

Phone: 1-800-363-1992 :

www.edmundoptics.com